
















AcuPulse 40WG Accessories

Part #	Description	Picture
FiberLase™ CO₂ Fibers		
AC-1059590	Single-use disposable, 2.0m long, 1.04mm OD, .5mm ID. Spot size: 320µm at fiber output. Up to 40 Watt input.	
Rigid Handpieces (silver handle)		
AC-1063930	FiberLase HP-R 60mm, straight, straight tip	
AC-1063970	FiberLase HP-R 180mm, straight, straight tip	
AC-1063920	FiberLase HP-R 60mm, straight, curved tip	
AC-1063940	FiberLase HP-R 140mm, straight, curved tip	
AC-1063960	FiberLase HP-R 180mm, straight, curved tip	
AC-1063980	FiberLase HP-R 240mm, bent, curved tip	
AC-1083980	FiberLase HP-R 140mm, bent, straight tip	
AC-1083970	FiberLase HP-R 240mm, bent, straight tip	
KT-1061370	FiberLase Rigid Handpiece Kit Includes all 8 rigid handpieces & FiberLase Handpiece Cleaning Kit (not shown)	(see individual images above)



Part #	Description	Picture
Malleable Handpieces (blue handle)		
AC-1083940	FiberLase HP-M 140mm, malleable	
AC-1083950	FiberLase HP-M 240mm, malleable	
AC-1083960	FiberLase Malleable Handpiece Kit: Includes both malleable handpieces & FiberLase Handpiece Cleaning Kit (not shown)	(see individual images above)
Accessories		
AC-1085400	FiberLase Handpiece Bending Tool	
SA-1065120	FiberLase Fiber Stripping Tool	
SA-1065040	FiberLase Fiber Cleaving Tool	
KT-1088480	FiberLase Fiber Renewal Kit: Includes Stripping, Cleaving & Magnification tool	(see individual images above)
SA-1065490	Sterilization Tray	
MD1475000	Internal Purge Air hose w/ Bacteriologic Filter Quantity: 10	