



VersaPulse® P20 Cart & Fiber Support Arm

Mobility & Convenience

The VersaPulse P20 holmium laser is a versatile, compact holmium laser that fits into a video cart. It provides the optimal power for effective stone and soft tissue treatments in a portable design that is perfect for a physician's office or small operating room.



VERSAPULSE P20 CART

The P20 cart provides the perfect mobile platform for your P20 and includes a convenient storage drawer.

- Durable aluminum frame offers high level of stability
- Drawer for storing fibers, goggles and other accessories
- Ergonomic design with push handles and locking wheels

VERSAPULSE P20 FIBER SUPPORT ARM

The VersaPulse P20 fiber arm helps ensure that your fiber is kept in the optimum position.

- Eliminates undue strain on fiber connection
- Keeps fiber elevated above the floor
- Adjustable joints and telescopic arm provide optimal support for the fiber



VersaPulse P20 cart



VersaPulse P20 fiber support arm



Lumenis Global Sales Offices

America

Santa Clara, CA

Tel + 1 408 764 3000
 + 1 877 586 3647
 Fax + 1 408 764 3999

Europe, Middle East & Africa

**Dreieich-Dreieichenhain,
 Germany**

Tel + 49 6103 8335 0
 Fax + 49 6103 8335 300

Formello (RM), Italy

Tel + 39 06 90 75 230
 Fax + 39 06 90 75 269

London, UK

Tel + 44 20 8324 4200
 Fax + 44 20 8324 4222

Asia-Pacific

Beijing, China

Tel + 86 10 6510 2620
 Fax + 86 10 6510 2621

Hong Kong

Tel + 852 2174 2800
 Fax + 852 2722 5151

Gurgaon, India

Tel + 91 124 422 33 80
 Fax + 91 124 422 33 83

Japan

Tokyo, Japan

Tel + 81 3 5789 8300
 Fax + 81 3 5789 8301

Internet: www.lumenis.com

Email: information@lumenis.com

Specifications

The VersaPulse P20 cart is fitted with 4 twin wheel swivel castors. The front 2 castors have locks. The cart also come with a lower shelf, push handles and a handy drawer for keeping your accessories.

VersaPulse P20 Cart Dimensions	
Part Number	SA-1069670
Width	714 mm
Depth	625 mm
Height	829 mm
Drawer Width	600 mm
Drawer Depth	575 mm
Drawer Height	96 mm
Load Capacity	150 kg / 330 lbs
Finish	RAL 7035 light grey powder-coated w/blue

Fiber Arm Support Dimensions	
Folded Length	235 mm
Extended Length	1230 mm



VersaPulse P20 cart drawer detail



VersaPulse P20 fiber support arm folded

