


SlimLine™ Fibers



Choose the size fiber wanted for any procedure

The family of SlimLine end-firing fibers is available in 200, 365, 550, and 1000 micron sizes for the broadest range of clinical applications. Reusable, robust design lowers the cost per case.

- Blue jacket enhances visibility during endoscopic surgery
- Delivers holmium and Nd:YAG wavelengths through a single fiber for convenience and utility
- Compatible with cystoscopes and ureteroscopes, including the smallest 7.5F flexible design



Fibers fit into the smallest flexible scope



Compatible with a wide range of cystoscopes and ureteroscopes

SlimLine™ Fibers

FORWARD FIRING FIBERS AND ACCESSORIES

SlimLine Fiber	200 micron	365 micron	550 micron	1000 micron
Product Number	0626-221-20	0624-148-36	0624-149-55	0624-497-01
Wavelength Compatibility	Holmium 45W Nd:YAG 45W	Holmium 100W Nd:YAG 100W	Holmium 100W Nd:YAG 100W	Holmium 100W Nd:YAG 100W
Scope Compatibility	Flexible, semi-rigid, & rigid endoscopes	Flexible, semi-rigid, & rigid endoscopes	Flexible, semi-rigid, & rigid endoscopes	Cystoscope & resectoscopes
Core	272 microns	365 microns	550 microns	940 microns
Maximum Outer Diameter	450 microns (1.35F)	580 microns (1.75F)	780 microns (2.25F)	1450 microns (4.5F)
Working Channel Needed	0.6 mm (1.8F)	.76 mm (2.3F)	1.0 mm (3.0F)	1.6 mm (5.0F)
Accessory Products				
Ceramic Scissors	5402-0152	5402-0152	5402-0152	5402-0152
Cleaving Tool	5402-0141	5402-0141	5402-0141	not used
Inspection Scope	0622-921-01	0622-921-01	0622-921-01	0622-921-01
Stripper	5402-0164	5402-0163	5402-0151A	5402-0159

FULL LINE OF ACCESSORIES

Accessories are available for stripping, cleaving, and sterilizing the reusable SlimLine Fibers.



VersaPulse® PowerSuite™ Holmium and Dual Wavelength Lasers



Sterilization Tray
Designed for processing fibers



Inspection Scope
Allows convenient visualization of the fiber tip and the connector



Fiber Stripper
Removes outer jacket from distal end of fibers

www.Lumenis.com

Lumenis Inc. - USA
Tel 800.635.1313
Fax 800.505.1133

Lumenis Ltd. - ISRAEL
Tel + 972.4.959.9000
Fax + 972.4.959.9050

Lumenis - HOLLAND
(BV Rep Office)
Tel + 010.65102620
Fax + 010.65102621

Lumenis SARL - FRANCE
Tel + 33.1.41.44.99.99
Fax + 33.1.47.28.99.09

Lumenis GmbH - GERMANY
Tel + 49.6071.968.110
Fax + 49.6071.968.400

Lumenis Ltd. - HONG KONG
Tel + 852.2174.2800
Fax + 852.2722.5151

Lumenis SARL - ITALY
Tel + 39.06.90.75.230
Fax + 39.06.90.75.269

Lumenis Co. Ltd. - JAPAN
Tel + 81.3.5805.5331
Fax + 81.3.5800.5730

Lumenis AB - SWEDEN
Tel + 46.44.21.21.26
Fax + 46.44.21.21.79

Lumenis Ltd. - UK
Tel + 44.20.8324.4200
Fax + 44.20.8324.4222

